



# TissueSpec® Bone ECM Hydrogel


Catalog # MTSBN101

TissueSpec® Bone ECM Hydrogel is a versatile extracellular matrix product comprised of bone-specific collagens and other ECM molecules of porcine origin. TissueSpec® hydrogels provide cells (e.g., osteoblasts, cancer cells, stem cells) a physiologic substrate for 3D cell culture that is easy to use and enhances cell function and cell-cell interactions.

## Features



- Derived from porcine bone tissue
- Contains bone-specific ECM components
- Supports primary cell and organoid cultures
- Compatible with standard cell culture protocols
- Consistent across lots
- Easy to use

## Applications in 3D cell culture




TissueSpec® Bone ECM Hydrogel can be applied as a thin gel to culture cells and study cellular activity. Cells interact with bone ECM and show bone-specific function. To study cell-cell interactions or microtissue structures, encapsulate cells or organoids within 3D TissueSpec® Bone ECM Hydrogel.

## TissueSpec® Bone ECM Hydrogel supports osteoblast phenotype



Primary human osteoblasts cultured in TissueSpec® Bone ECM Hydrogel (6 mg/mL) show **characteristic dendritic morphology** after 1 hour in culture. Osteoblasts have **higher alkaline phosphatase activity** in TissueSpec® Bone ECM Hydrogel than in Matrigel and collagen I after 7 days in culture.


# TissueSpec® Bone ECM Hydrogel enables cancer drug testing



TissueSpec® Bone ECM Hydrogel provides a model of breast cancer metastasis to bone, and shows **differential drug responses** of breast cancer subtypes compared to plastic (no ECM). Response to drug treatment (Paclitaxel, 5  $\mu$ M) and vehicle (DMSO) of (a) BT-549 cells and (b) T47-D cells for 48 hours. OD, optical density.

# TissueSpec® Bone ECM Hydrogel supports osteogenic differentiation

Primary human osteoblasts cultured in TissueSpec® Bone ECM Hydrogel support significant **calcium deposition (mineralization)** after 21 days in mineralization medium compared to osteoblasts cultured in standard growth medium and plastic (no ECM). Calcium deposition is confirmed by Alizarin red staining.



# TissueSpec® Bone ECM Hydrogel characteristics

## a Mass spec profile\*


ECM components	Biomolecules
collagens	type I, II, III, V, VI, X, XI, XII, XVI
glycoproteins	biglycan, osteonectin, periostin, tenascin C
proteoglycans	asporin, decorin, fibromodulin, lumican, osteoglycin, osteomodulin
matrix associated	albumin, annexin A2

\* partial list of components

## b Key components ( $\mu$ g/mL)


collagens (soluble)	7,500-9,500
elastin	400-500
glycosaminoglycans	100-200

## c Collagen



## d Mechanical properties

Bone ECM Hydrogel (6 mg/mL)	
Stiffness	55 $\pm$ 5 Pa



(a) Proteomic profile by mass spectrometry indicates that TissueSpec® Bone ECM Hydrogel has a unique, **bone-specific composition**. (b,c) TissueSpec® Bone ECM Hydrogel has a **consistent protein profile** across multiple lots. (d) Mechanical stiffness of TissueSpec® Bone ECM Hydrogel increases with hydrogel concentration and changes as cells remodel the bone matrix.